


Product Datasheet

Human TNFSF15 ELISA Kit (orb1293097)

Description	Human TNFSF15 (Tumor necrosis factor ligand superfamily member 15) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	TNFSF15
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O95150
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
TNFSF15(Tumor
necrosis factor
liga...