


Product Datasheet

Human TNFRSF18 ELISA Kit (orb1293098)

Description	Human TNFRSF18 (Tumor necrosis factor receptor superfamily member 18) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	TNFRSF18
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9Y5U5
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
TNFRSF18(Tumor
necrosis factor
rec...