


Product Datasheet

Human TNFRSF21 ELISA Kit (orb1293099)

Description	Human TNFRSF21 (Tumor necrosis factor receptor superfamily member 21) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	TNFRSF21
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	075509
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
TNFRSF21(Tumor
necrosis factor
rec...