


## Product Datasheet

### Human TYMS ELISA Kit (orb1293128)

|                        |                                             |
|------------------------|---------------------------------------------|
| <b>Description</b>     | Human TYMS (Thymidylate synthase) ELISA Kit |
| <b>Reactivity</b>      | Human                                       |
| <b>Range</b>           | 0.156-10ng/ml                               |
| <b>Target</b>          | TYMS                                        |
| <b>Storage</b>         | 2-8 °C for 6 months                         |
| <b>Note</b>            | For research use only                       |
| <b>Uniprot ID</b>      | <b>P04818</b>                               |
| <b>Sensitivity</b>     | 0.094ng/ml                                  |
| <b>Expiration Date</b> | Please enquire.                             |



Human  
TYMS(Thymidylate  
synthase) ELISA  
K...