


Product Datasheet

Human PIEZO1 ELISA Kit (orb1293139)

Description	Human PIEZO1 (Piezo-type mechanosensitive ion channel component 1) ELISA Kit
Reactivity	Human
Range	0.313-20ng/ml
Target	PIEZO1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q92508
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human
PIEZO1(Piezo-type
mechanosensitive...