


Product Datasheet

Human PNPLA2 ELISA Kit (orb1293142)

Description	Human PNPLA2 (Patatin-like phospholipase domain-containing protein 2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	PNPLA2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q96AD5
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human PNPLA2
(Patatin-like
phospholipase...)