


Product Datasheet

Human CLEC12A ELISA Kit (orb1293167)

Description	Human CLEC12A (C-type lectin domain family 12 member A) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	CLEC12A
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q5QGZ9
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human CLEC12A (C-type lectin domain family 12 member A) ELISA Kit