

## Product Datasheet

### Human CH25H ELISA Kit (orb1293171)

**Description** Human CH25H (Cholesterol 25-hydroxylase) ELISA Kit

**Reactivity** Human

**Range** 0.156-10ng/ml

**Target** CH25H


**Storage** 2-8 °C for 6 months

**Note** For research use only

**Uniprot ID** **O95992**

**Sensitivity** 0.094ng/ml

**Expiration Date** Please enquire.



Human  
CH25H(Cholesterol  
25-hydroxylase) ...