


Product Datasheet

Human AMDHD2 ELISA Kit (orb1293177)

Description	Human AMDHD2 (N-acetylglucosamine-6-phosphate deacetylase) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	AMDHD2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9Y303
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human AMDHD2
(N-
acetylglucosamine-
6-phos...