


Product Datasheet

Human ZBTB48 ELISA Kit (orb1293189)

Description	Human ZBTB48 (Zinc finger and BTB domain-containing protein 48) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	ZBTB48
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P10074
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
ZBTB48(Zinc
finger and
BTB domain-
...