


Product Datasheet

Human XRCC5/Ku80 ELISA Kit (orb1293191)

Description	Human XRCC5/Ku80 (X-ray repair cross-complementing protein 5) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	XRCC5/Ku80
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P13010
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
XRCC5/Ku80(X-
ray repair
cross-comp...