


Product Datasheet

Human UNC13A ELISA Kit (orb1293203)

Description	Human UNC13A/Munc13-1 (Protein unc-13 homolog A) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	UNC13A/Munc13-1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9UPW8
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
UNC13A/Munc13-
1(Protein unc-13
hom...