


Product Datasheet

Human TRIAP1 ELISA Kit (orb1293213)

Description	Human TRIAP1 (TP53-regulated inhibitor of apoptosis 1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	TRIAP1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O43715
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
TRIAP1(TP53-
regulated
inhibitor of...