


Product Datasheet

Human PSMA ELISA Kit (orb1293273)

Description	Human PSMA/FOLH1 (Glutamate carboxypeptidase 2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	PSMA/FOLH1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q04609
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
PSMA/FOLH1(Glutamate
carboxypeptid...