


Product Datasheet

Human PKNOX2 ELISA Kit (orb1293285)

Description	Human PKNOX2 (Homeobox protein PKNOX2) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	PKNOX2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q96KN3
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
PKNOX2(Homeobox
protein PKNOX2)
EL...