


Product Datasheet

Human PDE9A ELISA Kit (orb1293293)

Description	Human PDE9A (High affinity cGMP-specific 3',5'-cyclic phosphodiesterase 9A) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	PDE9A
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O76083
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
PDE9A(High
affinity
cGMP-
specific ...