

Product Datasheet

Guinea pig OVA sIgE ELISA Kit (orb1293320)

Description Guinea pig OVA sIgE (Ovalbumin Specific IgE) ELISA Kit

Reactivity Guinea pig

Range 1.563-100ng/ml

Target OVA sIgE


Storage 2-8 °C for 6 months

Note For research use only

Application notes Indirect ELISA Sandwich

Sensitivity 0.938ng/ml

Expiration Date Please enquire.



Guinea pig
OVA sIgE
(Ovalbumin
Specific ...