

Product Datasheet

Estrogen Receptor beta 2 peptide (orb12934)

Description	Estrogen Receptor beta 2 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	Estrogen Receptor beta 2
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb100531 (Estrogen Receptor beta 2 antibody)
Expiration Date	6 months from date of receipt.