

Product Datasheet

Mouse IgG1 kappa Isotype Control antibody (FITC) (orb1294012)

Description

Mouse monoclonal against Mouse

conjugated to FITC

Species/Host

Mouse

Conjugation

FITC

Tested

FC, IF

Applications
Preservatives

Mouse IgG1 kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.

Form/Appearance

Mouse IgG1 kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.

Storage

Shipped and store at 4°C for one year. Do not freeze.

Note

For research use only

Application notes

1. Take 100 µl peripheral blood anticoagulated by EDTA and add to the bottom of 5 ml tube. 2. Add 10 µl labeled antibody to the bottom of flow tube mixing with the whole blood, incubate for 20 minutes at room temperature away from light. 3. Add 2 ml RBC lysis buffer, incubate for 10 minutes away from light after mixing, dissolve red blood cells. 4. Sample tube is set to 1000 rpm centrifugation for 5 minutes, discard the supernatant. 5. Add 2 ml PBS wash buffer to resuspend the cells, then 1000 rpm centrifugation for 5 minutes, discard the supernatant. 6. Add 0.5 ml PBS wash buffer to resuspend the cells and detect by flow cytometry (sample should be determined on the day on the machine and can also be added fixation overnight at 4 °C then measured).

Clonality

Monoclonal

Source

Mouse

Dilution Range

10 µl / assay

Expiration Date

12 months from date of receipt.