



Product Datasheet

PPAR γ antibody (orb1294232)


Description	Rabbit polyclonal antibody to PPAR γ
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide
Target	PPAR γ
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500 – 1:1500 Immunofluorescence 1:300 – 1:
Expiration Date	12 months from date of receipt.



Anti-PPAR γ antibody at 1/1000 dilution,...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...