



Product Datasheet

CD81 Antibody (orb1294247)


Description	Rabbit polyclonal antibody to CD81
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / Extracellular domain corresponding to Human CD81
Target	CD81
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500 Immunohistochemistry (Paraffin) 1:100 - 1:200 Flow Cytometry 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Anti-CD81 antibody
at 1/500 dilution, Ly...



Immunohistochemical
analysis of paraffin...



Anti-CD81 antibody
at 1/500 dilution, Ly...