




Product Datasheet

Fibronectin 1/FN1 (C-term) Antibody (orb1294256)

Description	Rabbit polyclonal antibody to Fibronectin (C-term)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / encompassing a sequence within the C-terminus region.
Target	Fibronectin 1 / FN1 (C-term)
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300.
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500 Immunofluorescence 1:200-1:500 Immunohistochemistry (Paraffin) 1:100-1:200
Expiration Date	12 months from date of receipt.



All lanes: Anti-Fibronectin antibody at ...

Immunohistochemical
analysis of paraffin...Immunohistochemical
analysis of paraffin...