



Product Datasheet

Influenza A Virus Nucleoprotein Antibody (orb1294271)

Description	Rabbit polyclonal antibody to Influenza A Virus Nucleoprotein
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fragment corresponding to Influenza A Virus Nucleoprotein
Target	Influenza A Virus Nucleoprotein
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:6000-1:12000 Immunofluorescence 1:500 - 1:2000
Expiration Date	12 months from date of receipt.



All lanes:
Anti-
Influenza A
Virus
Nucleo...



Specificity
and
sensitivity
of Anti Infl...