

Product Datasheet

SARS-CoV-2 (COVID-19) nucleocapsid Antibody (orb1294283)

Description Rabbit polyclonal antibody to SARS-CoV 2 nucleocapsid protein

Species/Host Rabbit

Reactivity Virus

Conjugation Unconjugated

Tested Applications IF, NeA, WB

Immunogen Recombinant fragment corresponding to SARS-CoV 2 nucleocapsid protein

Target SARS-CoV-2 (COVID-19) nucleocapsid

Storage Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative


Note For research use only

Isotype Rabbit-IgG

Clonality Polyclonal


Dilution Range Western Blot 1:3000-1:5000
Immunofluorescence 1:500-1:1000

Expiration Date 12 months from date of receipt.



All lanes: Anti-SARS-CoV-2 nucleocapsid

...



Immunofluorescence:

Cells were stained

W...