



Product Datasheet

Toll-like Receptor 4/TLR4 Antibody (orb1294286)


Description	Rabbit polyclonal antibody to Toll-like Receptor 4
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / Intracellular domain corresponding to Human Toll-like Receptor 4 / TLR4
Target	Toll-like Receptor 4 / TLR4
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500–1:2000 Immunofluorescence 1:300–1:500 Immunohistochemistry (Paraffin) 1:200–1:300
Expiration Date	12 months from date of receipt.



All lanes: Anti-Toll-like Receptor 4 ant...



All lanes: Anti-Toll-like Receptor 4 ant...



Immunofluorescence: cells were fixed wit...