



## Product Datasheet

### Alix Antibody (orb1294288)


|                            |                                                                                                                                                                        |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ALIX                                                                                                                                     |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                 |
| <b>Reactivity</b>          | Human, Mouse, Rabbit                                                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                           |
| <b>Tested Applications</b> | IF, IHC, WB                                                                                                                                                            |
| <b>Immunogen</b>           | Synthetic peptide / encompassing a sequence within the center region                                                                                                   |
| <b>Target</b>              | Alix                                                                                                                                                                   |
| <b>Storage</b>             | Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative |
| <b>Note</b>                | For research use only                                                                                                                                                  |
| <b>Isotype</b>             | Rabbit-IgG                                                                                                                                                             |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                             |
| <b>Dilution Range</b>      | Western Blot 1:1000-1:2000 Immunofluorescence 1:400-1:600 Immunohistochemistry (Paraffin) 1:200-1:300                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                        |



All lanes: Anti-ALIX antibody at 1/1000 ...



Immunofluorescence: cells were fixed wit...



Immunohistochemical analysis of paraffin...