



Product Datasheet

Syndecan-1 antibody (orb1294301)


Description	Rabbit polyclonal antibody to Syndecan-1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / Extracellular domain corresponding to Human CD138 / Syndecan-1
Target	Syndecan-1 / CD138 Extracellular domain
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500-1:1000 Immunofluorescence 1:200-1:400 Immunohistochemistry (Paraffin) 1:200 - 1:300
Expiration Date	12 months from date of receipt.



All lanes:
Anti-
Syndecan-
1/CD138
antibod...



All lanes:
Anti-
Syndecan-
1/CD138
antibod...



All lanes:
Anti-
Syndecan-
1/CD138
antibod...