


Product Datasheet


Cytokeratin 7, CK7 Antibody (orb1294307)

Description	Rabbit polyclonal antibody to Cytokeratin 7
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / encompassing a sequence within the center region
Target	Cytokeratin 7, CK7
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:1000-1:2000 Immunofluorescence 1:200-1:500 Immunohistochemistry (Paraffin) 1:200-1:500
Expiration Date	12 months from date of receipt.




Anti-Cytokeratin 7 antibody - at 1/2000

...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...