



Product Datasheet

Pan-Cytokeratin, Pan-CK Antibody (orb1294309)


Description	Rabbit polyclonal antibody to pan Cytokeratin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Detects endogenous levels of total keratin 5, 6, 7, 8, 14, 17, 18, 19 and 20. The antibody does not cross-react with other keratins.
Target	Pan-Cytokeratin, Pan-CK
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C PBS (pH 7.4), 40% Glycerol. 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:2000–1:5000 Immunofluorescence 1:400–1:800 Immunohistochemistry (Paraffin) 1:400–1:800 Flow Cytometry 1:300–1:500
Expiration Date	12 months from date of receipt.



Anti-pan
Cytokeratin
antibody - at
1/300...



Anti-pan
Cytokeratin
antibody - at
1/300...



Cells were fixed with 4% PFA for 10 min ...