

Product Datasheet

EpCAM (Extracellular domain) antibody (orb1294311)

Description Rabbit polyclonal antibody to EpCAM antibody (Extracellular domain)

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications IF, IHC, WB

Immunogen Synthetic peptide / Extracellular domain corresponding to Human EpCAM

Target EpCAM / CD326 (Extracellular domain)

Storage Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative


Note For research use only

Isotype Rabbit-IgG


Clonality Polyclonal

Dilution Range Western Blot 1:1000-1:2000
Immunofluorescence 1:200-1:400
Immunohistochemistry (Paraffin) 1:100-1:150
Flow Cytometry 1:50-1:200


Expiration Date 12 months from date of receipt.



Anti-EpCAM antibody
- at 1/1000 dilution...



Immunohistochemical
analysis of paraffin...



Immunohistochemical
analysis of paraffin...