

Product Datasheet

EpCAM (intracellular domain) antibody (orb1294312)

Description Rabbit polyclonal antibody to EpCAM antibody (intracellular domain)

Species/Host Rabbit

Reactivity Human, Mouse, Rabbit

Conjugation Unconjugated

Tested Applications IF, IHC, WB

Immunogen Synthetic peptide / Intracellular domain corresponding to Human EpCAM

Target EpCAM / CD326 (Intracellular domain)

Storage Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative


Note For research use only

Isotype Rabbit-IgG


Clonality Polyclonal

Dilution Range Western Blot 1:1000-1:2000
Immunofluorescence 1:200-1:500
Immunohistochemistry (Paraffin) 1:200-1:300


Expiration Date 12 months from date of receipt.



Anti-EpCAM antibody
- at 1/1000 dilution...



Immunohistochemical
analysis of paraffin...



Immunohistochemical
analysis of paraffin...