



Product Datasheet

p63 Antibody (orb1294319)


Description	Rabbit polyclonal antibody to p63
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / encompassing a sequence within the N-terminus region.
Target	p63
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 1%BSA, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500-1:1000 Immunofluorescence 1:100-1:300 Immunohistochemistry (Paraffin) 1:100
Expiration Date	12 months from date of receipt.



All lanes: Anti-p63 antibody - at 1/500 ...



Immunohistochemical analysis of paraffin...



Immunofluorescence: cells were fixed wit...