



Product Datasheet

PERK/EIF2AK3 Antibody (orb1294328)


Description	Rabbit polyclonal antibody to PERK
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / encompassing a sequence within the center region
Target	PERK / EIF2AK3
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500 Immunofluorescence 1:300-1:500 Immunohistochemistry (Paraffin) 1:100-1:200
Expiration Date	12 months from date of receipt.



All lanes:
Anti-PERK
antibody - at
1/500...



All lanes:
Anti-PERK
antibody - at
1/500...



All lanes:
Anti-PERK
antibody - at
1/500...