



Product Datasheet

VEGF Antibody (orb1294346)


Description	Rabbit polyclonal antibody to VEGF
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / encompassing a sequence within the center region
Target	VEGF
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500-1:1000 Immunofluorescence 1:200-1:500 Immunohistochemistry (Paraffin) 1:100-1:300 Immunoprecipitation : Reported in the literature, not verified in house.
Expiration Date	12 months from date of receipt.



VEGF antibody at 1/1000 dilution, Lysate...



VEGF antibody at 1/1000 dilution, Lysate...



Immunofluorescence: cells were fixed wit...