



Product Datasheet

Cytochrome C/CYCS Antibody (orb1294348)


Description	Rabbit polyclonal antibody to Cytochrome C
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / encompassing a sequence within the center region
Target	Cytochrome C / CYCS
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:1000-1:2000 Immunofluorescence 1:200-1:500 Immunohistochemistry (Paraffin) 1:200-1:400
Expiration Date	12 months from date of receipt.



All lanes: Anti-Cytochrome C antibody at...



Immunohistochemical analysis of paraffin...



Immunofluorescence: cells were fixed wit...