

Product Datasheet

Caspase-9 cleaved Asp315 Antibody (orb1294352)

Description Rabbit polyclonal antibody to Caspase 9 cleaved Asp315

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications IF, IHC, WB

Immunogen Synthetic peptide corresponding to residues surrounding Asp315 of human Caspase 9

Target Caspase-9 cleaved Asp315

Storage Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative


Note For research use only

Isotype Rabbit-IgG


Clonality Polyclonal

Dilution Range Western Blot 1:500 Immunofluorescence 1:300-1:400 Immunohistochemistry (Paraffin) 1:100-1:300

Expiration Date 12 months from date of receipt.




All lanes: Anti-Caspase 9 cleaved Asp315...



All lanes: Anti-Caspase 9 cleaved Asp315...

Caspase-9



Immunofluorescence: cells were fixed wit...