


Product Datasheet

Caspase 9 cleaved Asp330 Antibody (orb1294353)

Description	Rabbit polyclonal antibody to Caspase 9 cleaved Asp330
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide corresponding to residues surrounding Asp330 of human Caspase 9
Target	Caspase 9 cleaved Asp330
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500 Immunofluorescence 1:300-1:400
Expiration Date	12 months from date of receipt.



Anti-Caspase 9 cleaved Asp330 antibody -...