



Product Datasheet

CD45 (Intracellular domain) antibody (orb1294366)


Description	Rabbit polyclonal antibody to CD45 antibody (Intracellular domain)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Synthetic peptide / Intracellular domain corresponding to Human CD45
Target	CD45 (Intracellular)
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500-1:1000 Immunofluorescence 1:200-1:500 Immunohistochemistry (Paraffin) 1: 200-1:400 Flow Cytometry 1:50-1:100
Expiration Date	12 months from date of receipt.



All lanes: Anti-CD45 antibody at 1/1000 ...



Product: CD45 antibody; Catalog: orb1294...



Immunohistochemical analysis of paraffin...