



Product Datasheet

IL-6 Antibody (orb1294370)


Description	Rabbit polyclonal antibody to IL6
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fragment corresponding to Human IL-6.
Target	IL-6
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500-1:700 Immunofluorescence 1:300-1:400 Immunohistochemistry (Paraffin) 1:100-1:200
Expiration Date	12 months from date of receipt.



Anti-IL6 antibody -
at 1/500 dilution,
L...



Anti-IL6 antibody -
at 1/500 dilution,
L...



Cells were fixed
with 4%
paraformaldehyd...