



Product Datasheet

Tau (Phospho-Thr181) antibody (orb1294393)


Description	Rabbit polyclonal antibody to Tau (Phospho Thr181)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, WB
Immunogen	Synthetic peptide corresponding to residues surrounding Thr181 of human Tau protein
Target	Phospho-Tau (Thr181)
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 30% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500-1:1500 Immunofluorescence 1:200-1:400 Immunohistochemistry (Paraffin) 1:100
Expiration Date	12 months from date of receipt.



Anti-Phospho-Tau (Thr181) antibody at 1/...



Immunohistochemical analysis of paraffin...



Specificity and sensitivity of Phospho-T...