



Product Datasheet

Tau Antibody (orb1294394)


Description	Rabbit polyclonal antibody to Tau
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P, WB
Immunogen	Synthetic peptide / encompassing a sequence within the C-terminus region.
Target	Tau
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:3000 – 1:5000 Immunofluorescence 1:200 – 1:500 Immunohistochemistry (Paraffin) 1:200 – 1:300 Immunohistochemistry (Frozen) 1:200 – 1:400
Expiration Date	12 months from date of receipt.



Anti-Tau antibody,
Lysates at 40 µg
per...



Immunohistochemical
analysis of paraffin...



Immunofluorescent
analysis. Sample:
prim...