



Product Datasheet

CD133 antibody (orb1294396)


Description	Rabbit polyclonal antibody to CD133
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / Extracellular domain corresponding to Human CD133.
Target	CD133 / Prominin-1 (Intracellular domain)
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500-1:1000 Immunofluorescence 1:300-1:500 Immunohistochemistry (Paraffin) 1:200-1:400
Expiration Date	12 months from date of receipt.



All lanes:
Anti-CD133
antibody at
1/700 ...



All lanes:
Anti-CD133
antibody at
1/700 ...



All lanes:
Anti-CD133
antibody at
1/700 ...