



Product Datasheet

Neurofilament L antibody (orb1294400)


Description	Rabbit polyclonal antibody to Neurofilament L
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Recombinant fragment corresponding to Human Neurofilament L. Encompassing a sequence within the N-terminus region.
Target	Neurofilament L / NF L
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 1% BSA, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:1000-1:2000 Immunofluorescence 1:200-1:500 Immunohistochemistry (Paraffin) 1:50-1:100
Expiration Date	12 months from date of receipt.



Anti-Neurofilament-L/NF-L antibody at 1/...



Immunofluorescent analysis. Sample: prim...



Immunohistochemical analysis of paraffin...