



## Product Datasheet

### Beta Tubulin 3 antibody (orb1294407)


|                            |                                                                                                                                                                                |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Beta Tubulin 3                                                                                                                                   |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                         |
| <b>Reactivity</b>          | Human, Mouse, Rabbit                                                                                                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                   |
| <b>Tested Applications</b> | IF, IHC-Fr, IHC-P, WB                                                                                                                                                          |
| <b>Immunogen</b>           | Synthetic peptide / encompassing a sequence within the C-terminus region.                                                                                                      |
| <b>Target</b>              | Beta Tubulin 3 / Tuj1                                                                                                                                                          |
| <b>Storage</b>             | Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 1% BSA, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative |
| <b>Note</b>                | For research use only                                                                                                                                                          |
| <b>Isotype</b>             | Rabbit-IgG                                                                                                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                     |
| <b>Dilution Range</b>      | Western Blot 1:5000 - 1:10000<br>Immunofluorescence 1:200-1:500<br>Immunohistochemistry (Paraffin) 1:100-1:300<br>Immunohistochemistry (Frozen) 1:100-1:300                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                |



All lanes: Anti-Beta Tubulin 3/Tuj1 anti...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis (Frozen sec...