



Product Datasheet

Nestin (NES) antibody (orb1294409)


Description	Rabbit polyclonal antibody to Nestin antibody (NES)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide / encompassing a sequence within the N-terminus region.
Target	Nestin / NES
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 1% BSA, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:500-1:1000 Immunofluorescence 1:200-1:500 Immunohistochemistry (Paraffin) 1:100-1:200
Expiration Date	12 months from date of receipt.



Anti-Nestin antibody
at 1/1000 dilution,...



Sample: Neuron 2A;
Nestin antibody
(orb1...



Immunohistochemical
analysis of paraffin...