



Product Datasheet

Vinculin Antibody (orb1294416)


Description	Rabbit polyclonal antibody to Vinculin
Species/Host	Rabbit
Reactivity	Human, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic peptide / encompassing a sequence within the center region.
Target	Vinculin
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Isotype	Rabbit-IgG
Clonality	Polyclonal
Dilution Range	Western Blot 1:2000 - 1:4000 Immunofluorescence 1:200 - 1:500 Immunohistochemistry (Paraffin) 1:200 - 1:300
Expiration Date	12 months from date of receipt.



Immunofluorescence:
cells were fixed wit...



All lanes: Anti-
Vinculin antibody -
Load...



All lanes: Anti-
Vinculin antibody -
Load...