



Product Datasheet

AKT (phospho-T450) antibody (orb129508)

Description	Rabbit polyclonal antibody to AKT (phospho-T450)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T450 of human Akt1(NP_005154.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	60kDa
Uniprot ID	P31749
Dilution Range	WB: 1:500-1000, IF/ICC: 1:20-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of MCF7 cell...