



Product Datasheet

APEX1 antibody (orb129513)


Description	Rabbit polyclonal antibody to APEX1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human APEX1/APE1(NP_542380.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39kDa
Uniprot ID	P27695
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa, HepG2 cells...



IHC-P staining of human lung cancer tissue...



Immunohistochemical analysis of formalin-fixed tissue...