


# Product Datasheet

## Bip antibody (orb129538)

|                            |                                                      |
|----------------------------|------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to Bip.                            |
| <b>Reactivity</b>          | Human, Mouse, Rat                                    |
| <b>Conjugation</b>         | Unconjugated                                         |
| <b>Tested Applications</b> | IHC, WB                                              |
| <b>Immunogen</b>           | A synthetic peptide of human BiP/GRP78(NP_005338.1). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide,50% glycerol,pH7.3.      |
| <b>Concentration</b>       | 1 mg/ml                                              |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.          |
| <b>Note</b>                | For research use only                                |
| <b>Isotype</b>             | IgG                                                  |
| <b>Clonality</b>           | Polyclonal                                           |
| <b>MW</b>                  | 78KDa                                                |
| <b>Uniprot ID</b>          | <b>P11021</b>                                        |
| <b>Dilution Range</b>      | WB: 1:500-2000                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                      |



Western blot  
analysis of  
extracts of  
var...