



Product Datasheet

CD133 antibody (orb129549)


Description	Rabbit polyclonal antibody to CD133, which is also known as prominin or AC133. CD133 is a marker fre...
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	A synthetic peptide of human CD133(NP_006008.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120KDa
Uniprot ID	043490
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of SW620 cells.



Immunohistochemical staining of human prostate tissue.



Immunohistochemical staining of human lung tissue.